



Instructions for Use

PHILADELPHIA[®]

Tracheotomy Collar

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	5
FR	Notice d'utilisation	6
ES	Instrucciones para el uso	7
IT	Istruzioni per l'uso	8
NO	Bruksanvisning	9
DA	Brugsanvisning	10
SV	Bruksanvisning	11
FI	Käyttöohjeet	12
NL	Gebruiksaanwijzing	13
PT	Instruções de Utilização	14
PL	Instrukcja użytkownika	15
RU	Инструкция по использованию	16
JA	取扱説明書	17
ZH	中文说明书	18
KO	사용 설명서	20

Size Chart				
Measurement A - Neck Height				
Infant	1.25" (3.2 cm)			
Pediatric	1.75" (4.4 cm)			
Adult	2.25" (5.7 cm)	3.25" (8.3 cm)	4.25" (10.8 cm)	5.25" (13.3 cm)

Measurement B - Neck Circumference	
Size	Circumference
Infant	6"–8" (15.2 – 20.3 cm)
Pediatric	8"–11" (20.3 – 27.9 cm)
Adult Small	10"–13" (25.4 – 33.0 cm)
Adult Medium	13"–16" (33.0 – 40.6 cm)
Adult Large	16"–19" (40.6 – 48.3 cm)
Adult Extra Large	19"+ (48.3 cm+)

ENGLISH

MEASUREMENT & SIZING

Proper sizing is important for patient immobilization and comfort. Always maintain the patient's head in neutral alignment.

Height:

Measure the vertical distance from the top of the patient's shoulder to the tip of the chin (**Figure 1-a**). The measurement corresponds to the HEIGHT of the collar. If the patient's measurement falls between two consecutive sizes, apply the smaller size first.

Circumference:

Measure around the patient's neck (**Figure 2-b**). The measurement corresponds to the size CIRCUMFERENCE of the collar.

APPLICATION

Proper application of the collar is as important as proper sizing for immobilization and patient comfort. The two-piece design consists of a front piece and a back piece which are packaged as a set.

1. After selecting the properly sized collar, apply the back piece of the collar to the back of the patient's neck. Center the collar. The back arrow should point upward.
2. Apply the FRONT piece of the collar with the chin secured in the recess. Center the collar to secure neutral alignment. The front piece OVERLAPS the back piece to ensure effective immobilization and comfort. The front arrow should point upward.
3. With hook and loop fasteners, tighten the collar with a bilateral adjustment. This will secure the patient's cervical region in neutral alignment.

FEATURES

- Non-toxic Plastazote® Foam
- Total Cervical Arch Support
- Full range of sizes for custom fit
- X-Ray Lucent
- Easily cleaned with soap and water
- Hook and loop fastener
- Use as prescribed by a physician

ADDITIONAL INSTRUCTIONS

- If you experience sudden neck or back pain, numbness or tingling, contact your doctor.
- If you have any problems or questions about the use and care of your brace, contact your doctor or orthotic professional.
- These instructions do not supercede the specific orders of your physician.

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