

















Quickie Built-4-Me

Our Built-4-Me Team is dedicated to creating product customizations beyond what is listed on our product order forms. Experts in the engineering and manufacturing of custom rehab products will listen to your specific requests and build a chair as individual as your fingerprint!

JAY[®] Seating Selection

JAY CUSHIONS

				
Basic E2601, E2602	Soft Combi™ P E2605, E2606	Care™ E2607, E2608	Zip™ E2607	
				
GO™ E2603, E2604	Ion™ E2607, E2608	Easy™ E2607, E2608	Lite™ E2603, E2604, E2607, E2608	Union™ E2607, E2608
				
X2™ E2622	J2™ E2622, E2623, E2624, E2625	J2™ Deep Contour E2622, E2623, E2624, E2625	Fusion™ E2622, E2623	J3™ E2622, E2623, E2624, E2625

JAY BACKS

		
Basic E2611	GO™ E2611, E2612	Care™ E2615
		
J2™ E2613	J2™ Deep Contour E2620	Focus Point™ E2620
		
J3 Posterior E2613, E2614	J3 Posterior & Lateral E2615	J3 Posterior Deep E2620



LIVE WITHOUT LIMITS.



		RIGID					TILT-IN-SPACE			FOLDING							
		K0009	K0005				K0009 & E1161	E1161		K0005					K0004	K0009	
		Q7 Active Rigid	Q7 Adjustable	GT	QRi	GP / GPV	GPSA	IRIS / IRIS SE	SR45	TS / TS SE	Quickie 2 Lite	Quickie 2	Quickie 2 HP	QXi	LXI	LX	M6
Specifications subject to change without notice.																	
Design	Frame Style	Rigid - Open	Rigid - Open	Rigid - Open	Rigid - Open	Rigid - Traditional	Rigid - Swing-Away	Rotation in Space	Rotation in Space	Tilt	Folding	Folding	Folding	Folding	Folding	Folding	Folding
	Standard Weight†	17.4 lbs (Std. Config.) 9.7 lbs (Transport)	18.4 lbs (Std. Config.) 11.1 lbs (Transport)	19.6 lbs	22 lbs	23.75 lbs (GP) 23.6 lbs (GPV)	26 lbs	39 lbs (w/footrests, 18x18)	39 lbs	49 lbs	25.5 lbs (fixed) 27.75 (swing-away)	28.9 lbs	27.9 lbs	28.8 lbs	29.6 lbs	29.6 lbs	53 lbs (w/o footrests)
	User Weight Capacity	265 lbs	265 lbs	265 lbs	265 lbs	250 lbs*	250 lbs*	250 lbs (Std) 350 lbs (HD)	265 lbs	250 lbs (Std) 350 lbs (HD)	265 lbs	265 lbs - STD 350 lbs - HD	250 lbs	265 lbs.	265 lbs	250 lbs	650 lbs
Adjustability	Seat Width	12" to 20"	12" to 20"	12" to 20" *	12" to 20"	12" to 22" *	12" to 22" *	14" to 22"	14" to 20"	14" to 20"	14" to 20"	12" to 22" *	11" to 22"	14" - 20"	12" to 20"	12" to 20"	22" to 30"
	Seat Depth	12" to 20"	12" to 20"	12" to 22" *	12" to 20"	12" to 18" * (GP) 12" to 22" * (GPV)	12" to 22"	14" to 22"	14" to 22"	14" to 22"	14" to 20"	12" to 20" *	10" to 20"	14" - 20"	12" to 20" *	12" to 18"*	18" to 22"
	Back Angle Adjustment	Fixed	-8° to 27°	-8° to 17°	-8° to 17°	-4° to 8°	-4° to 8°	-30° to 5°	-30° to 5°	-30° to 12°	-14° to 4°	-15° to 3°	-15° to 3°	-14° to 5°	-14° to 5°	-14° to 5°	N/A
	Camber	0°, 3°, 6°	0°, 3°, 6°	0°, 3°, 6°, 9°	0° to 3°	0° to 11°	0° to 11°	0°	0°	0°, 4°	0°, 2°, 4°, 6°	0°, 2°, 4°, 6°	0° to 7°	0°, 3°	1° to 5°	1° to 5°	0°
	Tension Adj. Upholstery	Yes	Yes	Yes	Yes	Yes	Yes	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	No
Clinical Applications	Front Seat-to-Floor Ht.	16" to 21"	16" to 21"	16" to 21"	15" to 20"	15.5" to 23"	15.5" to 23"	12.5" to 19"	14" to 20.5"	—	15.5" to 21"	14" to 22"	17" to 22"	12.5" - 19.5"	16" to 20.5"	16" to 20.5"	17" to 20"
	Rear Seat-to-Floor Ht.	13" to 20"	13" to 20"	13" to 20"	14" to 20"	14.5" to 20.5" (GP) 14.5" to 20.6" (GPV)	14.5" to 20.5"	12.5" to 19"	14" to 20.5"	—	13" to 19.5"	13" to 22"	15.5" to 21"	12" - 19"	15" to 19.5"	15" to 19.5"	17" to 20"
	Footrest/Frame Angle	75°, 80°, 85°	70°, 80°, 85°	75°, 80°, 85°	75°, 80°, 85°	70°, 80°, 85°	Swing-Away	60°, 70°, 80°, 90°, Elevating	70°, 80°, Elevating	60°, 70°, 80°, 90°, Elevating	70°, 80°	Swing-Away	60°, 70°	60°, 70°, 80°, Elevating	Swing-Away	Swing-Away	Swing-Away
	Fold-Down Back	No	Yes	Yes	Yes	Yes	Yes	Yes	Optional	No	No	No	No	No	No	No	No
	Adjustable Back Height	Fixed	Adj. 8" to 20"	Adj. 10" to 20"	Adj. 14" to 20"	Adj. 8" to 20"	Adj. 8" to 20"	Adj. 15" to 21" or 18" to 24"	Adj. 15" to 21" or 18" to 24"	No	9" to 21"	9" to 21"	8.5" to 20"	13" - 17" or 16" - 20"	14" to 19"	14" to 19"	16.75" to 20.13"
	Center of Gravity (COG)	0.5" to 4"	1" to 4"	-1" to 4.5"	1" to 3"	0.5" to 3.5"	0.5" to 3.5"	N/A	N/A	N/A	0.5" to 3"	-1" to 4.5"	0.5" to 3.5"	1" to 3"	0.5" to 3.5"	0"	N/A
	Tilt Range	—	—	—	—	—	—	40° or 55°	45°	53°	—	—	—	—	—	—	—
Lifestyle Options	Xtender™ Power Assist	No	Yes	Yes	No	Yes	Yes	No	No	No	No*	Yes	Yes	No	No	No	No
	Front Suspension	Yes	Yes	Yes	Yes	Built-4-Me	Built-4-Me	Built-4-Me	No	No	Yes	Yes	Yes	No	Built-4-Me	Built-4-Me	No
	Rear Suspension	No	Yes	Yes	No	No	No	No	No	No	No	No	No	No	No	No	No
	Transit	Yes	Yes	No	Yes	No	No	Yes	Yes	No	Yes	Yes	Yes	No	Yes	No	No
	Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

* Additional options available through Quickie Built-4-Me. † Transport Weight = all items that can be removed without tools to depict transporting into or out of a vehicle. Standard Weight = 16" x 16" chair with wheel locks.

Quickie Colors



IMPORTANT NOTE: The color options shown are for reference only. They are not intended as accurate representations. Refer to the Quickie Color Brochure or a Quickie Color Ring for representations of Quickie colors.

IMPORTANT NOTE: Print processes cannot accurately represent paint finishes. Actual colors may vary.

*Black Zebra Kolorfusion™ pattern is available on any base color. Shown here on Candy Blue base color.