EIO Push Chair 2 year warranty




Affordable Chair and Positioning

- Designed for children from early intervention to pre-adolescent ages
- Easy to fold and transport
- Lightweight aluminum frame
- Calf panel and seat back lay your child flat
- Standard features: see next page

Consider Towards FIO Deads Chair	i. 6300000	
Special Tomato [®] EIO Push Chai Stroller Only Seat to Foot Well	13"	
Liners Size 1 and 2 Seat to Foot Well	15"	
Sitter Size 1 Seat to Foot Well	16"	
Sitter Size 2 Seat to Foot Well	17"	
Open Chair Dimensions	45"H × 24"W × 46"D	
Folded Chair Dimensions	16"H × 18"W × 32"D	95 Tob
Handle Height to Floor	30" - 43"	
Seat	14"W × 12"D	
Back Height	24"	
Max. User Weight	90 lbs.	
Stroller Weight	22 lbs.	
Includes foot-bloc	2-position ck.	



Features & Accessories



- 1 Size 1 and 2 Soft-Touch® Liners
 Attach to the stroller, with color choice
 of Cherry, Lilac, Aqua or Dark Gray.
 Sold separately, see page 16.
- 4 Shopping Basket 10 lb. weight capacity.
- Sun Shade Adjustable push handle, padded headrest and removable sun canopy with window, which includes a large pocket.

- Size 1 and 2 Soft-Touch® Sitters Attach to the stroller, with color choice of Lilac, Aqua or Dark Gray. Sold separately, see page 9.
- 5 Angle Adjustable Seat Back (25° to 75°) and foot well (+15° to -55°).
- 8 Rain Cover 63030000 Sold separately.



- Adjustable 8-point Harness
 5-Point comfort trunk harness with padded shoulders and chest strap.
 Also includes a padded 3-point pelvic harness for hip and lateral support.
- Solid Filled Wheels Solid filled 9" front and 11" rear wheels require no air. Swivel front wheels with lock-out feature add security on rough terrain. Includes one-touch parking brake.