

LEVO Kid & Junior

The every day wheelchair that let's
your child stand, play and grow –
everyday

LEVO Kid & Junior



LEVO stands for a better quality of life

LEVO.[®] **dt** 
THE EXPERTS IN STANDING

No limitations for a life full of promise!



To be part of...


Your child is active and wants to reach up and be part of it all. The LEVO Kid & Junior is their partner in play, in learning, in growing to their full potential. Like your child, the LEVO Kid & Junior is lively and fun with special abilities. Different seating positions and standing up to 86° offers a high degree of independence and flexibility. Easy adjustable growth plates, a wide range of accessories and options let the LEVO Kid & Junior smoothly grow with your child for many years.



LEVO Stands for a better quality of life

Highlights

- Smooth standing and positioning movement
- Compact, maneuverable, light
- Safe and secure stand-up operation elevating to any and all angles between sitting and 86° standing
- Switch activation and switch option for all levels of hand/finger function
- Adjustable for upper and lower leg length, foot position and back height
- Six ground stability points ensure solid and secure standing position
- User weight up to 120 kg
- Wide range of different options and accessories
- Multiple color option
- Certificates: TÜV, CERAH, HI, FDA/510(k), Crash Test



Standard configuration

Armrests folding upwards and height and angle adjustable, folding down backrest angle adjustable, battery and charger, chest belt and knee support

Options/accessories LCEV

Foam seat cushion, Anti-tip rollers, Tip aid, Backrest extension, Push handles, Push handle extension, Swing away chest support, Spinerge wheels, Camber, Spoke guard, V-Trak back, ... many more on request

Technical Data

Model	LCEV		Kid & Junior	
	S - M	L - XL	Kid	Junior
Seat width (SW)	35/38/40/42/45/48/50		30/35	
Total width (W)	53/56/58/60/63/66/68		49/54	
Total length (L)	84 - 90	93 - 104	69 - 72	77 - 83
Total height (incl. Backrest) (H)	80		55	66
Backrest upholstery height min.	31		28	
Seat upholstery height front	48/51/54/57		33	40/44
Seat upholstery depth	29 - 42	43 - 57	24 - 40	
Footrest	Height and angle adjustable			
Total weight (w/o detachable parts)	25 (19)	26 (20)	17 (12)	19 (14)
Rear wheels	22" / 24" / 26"		16"	20" / 22"
Batteries	2 Batteries à 12V/2,7 Ah			
Transport dimensions max. (TLxTWxTH)	98 x 60 x 45		62 x 45 x 43	

Discrepancies in the measurements may arise due to wheelchair modifications. All measurements in cm and kg.

Options/accessories Kid & Junior

Growth kit seat width 30 to 35 cm, Growth kit Kid to Junior Height adjustable stroller handle, Attendant control for stand-up operation, Anti-tip devices, Table tray, Seat cushion, Back rest extension 6 to 16 cm, V-Trak back, Spinerge wheels, ... many more on request

Colors

5 standard colors (blue, red, yellow, platinum, black) and 190 other color options

Distributor:

Manufactured by:

LEVO AG
Switzerland

Anglikerstrasse 20
CH-5610 Wohlen

Tel. +41 (0)56 618 44 11
Fax +41 (0)56 618 44 10

office@levo.ch
www.levo.ch

LEVO.dt
A Dane Technologies Company