

Hill-Rom.
1000 Bed



Introducing the Hill-Rom® 1000 Bed.

Value minded doesn't mean sacrificing innovation and style. The Hill-Rom 1000 Bed fits the need for both patient and caregiver - dedicated to safety, comfort and efficiency.

At a glance

- One-button Dining Chair® position for patient comfort and safety
- Patient pendant is stowable in any siderail
- Shearless Pivot® mechanism with auto contour means less need for caregivers to reposition patients
- Higher headboard makes moving the bed easier
- Easy-to-read Line-of-Site® angle indicators for degree of head section and Trend and Reverse Trend
- Instant CPR release lever includes a dampened descent, gently lowering the patient and reducing caregiver back strain
- Point-of-Care® controls save caregiver steps
- Four dual-locking casters keep the bed in place
- Bed exit system
- Bed and mattress are designed to address main entrapment areas
- Optional built-in scale eliminates the need to move patients to a separate scale

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

Standard Features Unique to Hill-Rom

Patient Safety

Your patients' safety and comfort are your primary concern. The Hill-Rom® 1000 Bed incorporates several unique features to help you protect and care for your patients.



One-button Dining Chair® feature

- One-button control moves patient in and out of chair position.
- The bed maintains a low, safe height while putting patient in a chair position.



Stowable patient pendant

For ease of use, patients have full access to bed controls regardless of rail position.



Unique Brake Bar at the Footend of the Bed

- The brake bar is easily accessible even when siderails are down.
- The brake bar is located away from patient egress points.



Shearless Pivot® feature

This unique feature can reduce the need to reposition patients with the touch of a button.

Caregiver Ergonomics

With features that help provide safer patient mobility and ease of use, the Hill-Rom 1000 Bed was designed with caregiver safety as a priority.



Non-moving headboard

- Headboard remains at ergonomic height during all bed articulations.
- Stationary headboard height means no more concerns about IV poles breaking the lights.

Line of Site® angle indicator

The indicator is on the siderail for easy visibility.

Caregiver Convenience

The Hill-Rom 1000 Bed controls are created to be easy-to-use and promote efficient patient care.



Point-of-Care® controls

Identical, easy-to-read caregiver controls, including scale, are located on both sides of the bed.

Battery back-up

- There is no need to find an outlet when repositioning the bed during patient transfer.
- Battery back-up automatically charges when the bed is plugged in.

Additional Hill-Rom® 1000 Bed Advantages

The Hill-Rom 1000 Bed meets FDA Guidance on Hospital Bed System Dimensions (HBSW) to reduce the potential for siderail entrapment.

Scale

- Built-in scale eliminates the need for sling scale procedures helping to improve patient safety.
- Located on the siderail for convenient, easy-to-read controls.
- Scale can weigh in any position.

Bed Exit System

- Setting the bed exit system, combined with falls prevention protocols, can help to reduce patient falls.
- Newer, more reliable technology means fewer false calls.
- This system is single-mode with both tone and volume controls.



Technical Specifications

Hill-Rom 1000 Bed

Feature	Dimension	Feature	Dimension
Frame width.....	39.74" (100.95 cm)	Low position floor to deck.....	15.75" (42.55 cm)
Frame length with head and footboard (with bumpers) normal.....	90.04" (228.71 cm)	High position floor to deck.....	32.5" (82.55 cm)
Sleep deck width.....	36.00" (91.44 cm)	Max head elevation.....	65 degrees
normal length.....	80.00" (203.20 cm)	Max Trend and reverse Trend.....	16 degrees
Bed clearance under frame.....	5.00" (12.70 cm)	Caster size.....	6" (15 cm)
Optional Scale System accuracy.....	+/- 2.2 lbs or 1.1% of patient weight, whichever is greater		
weigh in any position.....	yes		
scale repeatability.....	2% of patient weight		

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

France +33 (0)2 97 50 92 12
 United Kingdom +44 (0)1530 411000
 Deutschland 0800 8 63 55 37
 Nederland +31 (0)347 32 35 32
 Italia + 39 (0)2 950541
 Suisse/Schweiz +41 (0)21 706 21 30
 (deutschsprachig) 0800 8 63 55 37

Österreich 0800 8 63 55 37
 Ireland +353 (0)1 413 6005
 Iberia +34 (0)93 6856000
 Nordic Countries +46 (0)20 78 10 30
 Export Contact your local Hill-Rom distributor or contact your Area Manager via website or call +1 812 934 8173

www.hill-rom.com

Hill-Rom® reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom® makes is the express written warranty extended on the sale or rental of its products.
 © 2013 Hill-Rom® SARL, ALL RIGHTS RESERVED.