



### SC-3008-DL

The ultimate pneumatic compression device. The SC-3008-DL Sequential Circulator is as unique as the patient population it services. All patients have different needs so our purpose-built technology measures and then delivers pressure to each and every chamber in the garment to give the best possible treatment.

The SC-3008-DL lets the user program the pressure setting in any and all of the 8 chambers, always keeping greater pressure distal so as to avoid any reflux. Complete with a compliance meter to monitor patient use, the SC-3008-DL also has an adjustable countdown clock that shuts off the pump after treatment has ended.

### Key Features:

- 8 Chamber System
- Gradient (8 levels)
- Sequential
- Non-peristaltic
- Distal to proximal cycle
- Unilateral or bilateral use
- 60 second cycle time
- 0 - 120 mmHg pressure range
- Slim and light design
- Silent
- Easy handling
- Timer
- Compliance meter
- Individual pressure adjustment in every chamber



### Specifications:

Dimensions: 4.5" x 11.75" x 7.75"  
11 cm x 30 cm x 20 cm

Weight: 5.9 lb  
2,68 kg

Supply Item: E0652



### SC-2008-OC

### Key Features:

- 8 Chamber System
- Gradient (8 levels)
- Sequential
- Non-peristaltic
- Distal to proximal cycle
- Unilateral or bilateral use
- 60 second cycle time
- 20 - 80 mmHg pressure range
- Slim and light design
- Silent
- Easy handling

### Specifications:

Dimensions: 4.5" x 11.75" x 8.125"  
11 cm x 30 cm x 21 cm

Weight: 5.25 lb  
2,38 kg

Supply Item: E0651